

● Three-dimensional analysis by TOF-SIM

TNE0006

Overview

In Time-of-Flight Secondary Ion Mass Spectrometry (TOF-SIMS), a separate sputter etching mirror (O_2 , Cs) is provided apart from the primary ion mirror. Analysis in the depth direction is performed by repeated sputter etching and normal measurement. While analysis in the depth direction can also be performed by dynamic SIMS, AES and XPS (ESCA) etc, TOF-SIMS offer the following advantages over the other methods:

- Allows data collection without restrictions on the number of components. (when both positive and negative polarities are measured)
- Allows three-dimensional imaging, and the extraction of two-dimensional images for any given location.
- Allows the extraction of a profile for any given location in the area under analysis.

Table1. Features of different methods

| | TOF-SIMS | FE-AES | XPS | SIMS |
|---|-----------------|------------|------------|-----------------|
| Detection sensitivity | ○ | △ | △ | ◎ |
| Addition of components after measurement | ○ | × | × | × |
| Minimum size and maximum depth | Several μm | <1 μm | 50 μm | Tens of μm |
| Under actual level conditions when measuring in depth direction | ↓ | | ↓ | ↓ |

Example: Three-dimensional analysis of TFT element




Fig.1 Depth profile.




Fig.2 Two-dimensional slice images (Cr^+).




Fig.3 Three-dimensional images.

Created by SCAS: <https://www.scas.co.jp/technical-informations/technical-news/pdf/tn339.pdf>

Total solution partner for analyses, measurements, and evaluations.

SCASS
SCAS SINGAPORE PTE LTD

[Laboratory] 17 Sakra Road, Pulau Sakra, Singapore 627886
[Office] 1 Gateway Drive, #09-09, Westgate Tower, Singapore 608531
E-mail: scass@scass.com.sg WEB site: <https://www.scass.com.sg/>